



Recombinant Human Uncharacterized protein C6orf141 (C6orf141)

Product Code	CSB-MP725558HU
Abbreviation	C6orf141
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SZD1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNDPFARMET RGPQGAANPM DSSRSLGDLG PFPREVGRGA PLAPGARNPA TAGASRSQGG GHEDRTADRA LGPRAGEELD RESWREKVL FLLHPERWLG TRGDPAREEV AGAEDLPHAG GEDHGEEPNY PSVFQRQKRI SGRRVAPPRD AADPPKYVLV RVEDYQVTQE VLQTSWAKGR MTTTRTEEHFV TALTFRSSRE GQPGERWGPA ESRALQARTG ASRVHAAGR R VSPSPGTWLE EIKL
Source	Mammalian cell
Target Names	C6orf141
Protein Names	Recommended name: Uncharacterized protein C6orf141
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.