



# Recombinant Human Protein FAM102A (FAM102A)

<b>Product Code</b>	CSB-BP725894HU
<b>Abbreviation</b>	FAM102A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5T9C2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAFLMKKKKF KFQTTFTLEE LTAVPFVNGV LFCKVRLLDG GDFVSLSSRE EVQENCVRWR KRFTFVCKMS ANPATGLLDP CVFRVSVRKE LKGGKAYSKL GFADLNLAEF AGSGSTVRCC LLEGYDTKNT RQDNSILKVT IGMFLLSGDP CFKTPPSTAK SISIPGQDSS LQLTCKGGGT SSGGSSTNSL TGSRPPKARP TILSSGLPEE PDQNLSSPEE VFHSGHSRNS SYASQQSKIS GYSTEHSRSS SLSDLTHRRN TSTSSASGG LGMTVEGPEG SEREHRPPEK PPRPPRPLHL SDRSFRRKKD SVESHPTWVD DTRIDADAIV EKIVQSQDFT DGSNTEDSNL RLFVSRDGSA TLSGIQLATR VSSGVYEPVV IESH
<b>Source</b>	Baculovirus
<b>Target Names</b>	FAM102A
<b>Protein Names</b>	Recommended name: Protein FAM102A Alternative name(s): Early estrogen-induced gene 1 protein
<b>Expression Region</b>	1-384
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.