



Recombinant Human Protein FAM102A (FAM102A)

Product Code	CSB-MP725894HU
Abbreviation	FAM102A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T9C2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAFLMKKKKF KFQTTFTLEE LTAVPFVNGV LFCKVRLLDG GDFVSLSSRE EVQENCVRWR KRFTFVCKMS ANPATGLLDP CVFRVSVRKE LKGGKAYSKL GFADLNLAEF AGSGSTVRCC LLEGYDTKNT RQDNSILKVT IGMFLLSGDP CFKTPPSTAK SISIPGQDSS LQLTCKGGGT SSGGSSTNSL TGSRPPKARP TILSSGLPEE PDQNLSSPEE VFHSGHSRNS SYASQQSKIS GYSTEHSRSS SLSDLTHRRN TSTSSASGG LGMTVEGPEG SEREHRPPEK PPRPPRPLHL SDRSFRRKKD SVESHPTWVD DTRIDADAIV EKIVQSQDFT DGSNTEDSNL RLFVSRDGSA TLSGIQLATR VSSGVYEPVV IESH
Source	Mammalian cell
Target Names	FAM102A
Protein Names	Recommended name: Protein FAM102A Alternative name(s): Early estrogen-induced gene 1 protein
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.