



Recombinant Human Cytochrome P450 4A22 (CYP4A22)

Product Code	CSB-BP725898HU
Abbreviation	CYP4A22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TCH4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	VLSPSR RLGGVSGILQ VTSLILLLLL LIKAAQLYLH RQWLLKALQQ FPCPPSHWLF GHIQEFQHDQ ELQRIQERVK TFPSACPYWI WGGKVRVQLY DPDYMKVILG RSDPKSHGSY KFLAPRIGYG LLLLNGQTWF QHRRMLTPAF HNDILKPYVG LMADSVRVML DKWEELLGQD SPLEVFQHVS LMTLDTIMKS AFSHQGSIQV DRNSQSYIQA ISDLNSLVFC CMRNAFHEND TIYSLTSAGR WTHRACQLAH QHTDQVIQLR KAQLQKEGEL EKIKRKRHLD FLDILLAKM ENGSILSDKD LRAEVDTFMF EGHDTTASGI SWILYALATH PKHQERCREE IHGLLDGAS ITWNHLDQMP YTTMCIKEAL RLYPPVPGIG RELSTPVTFP DGRSLPKGIM VLLSIYGLHH NPKVWPNLEV FDPSRFAPGS AQHSHAFLPF SGGSRNCIGK QFAMNQLKVA RALTLLRFEL LPDPTRIPIP MARLVLKSKN GIHLRLRRLP NPCEDKDQL
Source	Baculovirus
Target Names	CYP4A22
Protein Names	Recommended name: Cytochrome P450 4A22 Alternative name(s): CYP1VA22 Fatty acid omega-hydroxylase Lauric acid omega-hydroxylase EC= 1.14.15.3
Expression Region	5-519
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This gene is part of a cluster of cytochrome P450 genes on chromosome 1p33.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.