



# Recombinant Human Uncharacterized protein C1orf114 (C1orf114)

<b>Product Code</b>	CSB-EP725903HU-B
<b>Abbreviation</b>	C1orf114
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5TID7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNENKDTDSK KSEEYEDDFE KDLEWLINEN EKSDASIIEM ACEKEENINQ DLKENETVME HTKRHSDPK SLQDEVSPRR NDIISVPGIQ PLDPISDS ENSFQESKLE SQKDLEEEED EEVRRYIMEK IVQANKLLQN QEPVNDKRER KLKFKDQLVD LEVPPLEDTT TFKNYFENER NMFGLKSLQC ISNDFGQEDV LLSLTNGSCE ENKDRILVE RDGKFELLNL QDIASQGFLP PINNANSTEN DPQQLPRSS NSSVSGTKKE DSTAKIHAVT HSSTGEPLAY IAQPPLNRKT CPSSAVNSDR SKGNGKSNHR TQSAHISPVT STYCLSPRQK ELQKQLEEKR EKLKREEERR KIEEKEKKR ENDIVFKAWL QKKREQVLEM RRIQRAKEIE DMNSRQENRD PQQAFRLWLK KKHEEQMKER QTEELRKQEE CLFFLKGTEG RERAFKQWLR RKRMEKMAEQ QAVRERTRQL RLEAKRSKQL QHHLYMSEAK PFRFTDHYN
<b>Source</b>	E.coli
<b>Target Names</b>	CCDC181
<b>Protein Names</b>	Recommended name: Uncharacterized protein C1orf114
<b>Expression Region</b>	1-509
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.