



Recombinant Human Uncharacterized protein C1orf114 (C1orf114)

Product Code	CSB-MP725903HU
Abbreviation	C1orf114
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TID7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNENKDTDSK KSEEYEDDFE KDLEWLINEN EKSDASIIEM ACEKEENINQ DLKENETVME HTKRHSDPK SLQDEVSPRR NDIISVPGIQ PLDPISDS ENSFQESKLE SQKDLEEEED EEVRRYIMEK IVQANKLLQN QEPVNDKRER KLKFKDQLVD LEVPPLEDTT TFKNYFENER NMFGLKSLQC ISNDFGQEDV LLSLTNGSCE ENKDRILVE RDGKFELLNL QDIASQGFLP PINNANSTEN DPQQLPRSS NSSVSGTKKE DSTAKIHAVT HSSTGEPLAY IAQPPLNRKT CPSSAVNSDR SKGNGKSNHR TQSAHISPVT STYCLSPRQK ELQKQLEEKR EKLKREEERR KIEEKEKKR ENDIVFKAWL QKKREQVLEM RRIQRAKEIE DMNSRQENRD PQQAFRLWLK KKHEEQMKER QTEELRKQEE CLFFLKGTEG RERAFKQWLR RKRMEKMAEQ QAVRERTRQL RLEAKRSKQL QHHLYMSEAK PFRFTDHYN
Source	Mammalian cell
Target Names	CCDC181
Protein Names	Recommended name: Uncharacterized protein C1orf114
Expression Region	1-509
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.