



Recombinant *Macaca fascicularis* Homeobox protein NANOG (NANOG)

Product Code	CSB-YP725910MOV
Abbreviation	NANOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TM84
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	>85% (SDS-PAGE)
Sequence	MSVDPACPQS LPCLEASDSK ESSPMPVICG PEENYPQLQM SSAEMPHTET VSPLPSSMDL LIQDSPDSST SPKGGKQPTAA ENSATKKEDK VPVKKQKART VFSSAQLCVL NDRFQRQKYL SLQQMQELSN ILNLSYKQVK TWFQNQRMKS KRWQKNNWPK NSNGVTQKAS APTYPSLYSS CHQGCLVNPT GNLPMWSNQT WNNSSWSNQT QNIQSWSNHS WNAQTWCTQS WNNQAWNSPF SNCGEESLQS CLQFQPNSPA SDLEAALEAA GEGLNVIQQT TRYLSTPQTV DLLLNYSTNM QPEDV
Source	Yeast
Target Names	NANOG
Protein Names	Recommended name: Homeobox protein NANOG Alternative name(s): Homeobox transcription factor Nanog
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.