



Recombinant Danio rerio Radial spoke head protein 9 homolog (rsph9)

Product Code	CSB-MP725914DIL
Abbreviation	rsph9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TYW6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDSDSLHYSL DLAAANGLTL SSEQRAALQS SLLIVKRNYK FSRVLFWGKI LGIKSDYFIA QGVEDDELKN RKSLSLNCV DWHLLPPATE SMIADVALAA TGRFTGDP SH EYEHEIRTE GEGDEATHEE VTVKVIASR LAIVSNIDK DVSVPARGAF IKSPNGKVQT NRSFGGLHPT EAAKLRNYLH FREPVNLRNK SILEMSELNP ALDFLDPLSE DILKGSWSLQ LDRGGTVCVL RLLWLGFTF FHVPQTPQHG YIYMGDGLMN LDLPFML
Source	Mammalian cell
Target Names	rsph9
Protein Names	Recommended name: Radial spoke head protein 9 homolog
Expression Region	1-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.