



Recombinant Danio rerio BTB/POZ domain-containing protein 6-B (btbd6b)

Product Code	CSB-BP725916DIL
Abbreviation	btbd6b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TZE1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAAELYPASI NTNLPNSNST AVTAASKKTI VQVTQTVTTP TTTATQQNIN NNNVETASWQ STHPTLRERN ALMFNNELMA DVHFVVGPPG ASQKVPAHKY VLAVGSSVFG AMFYGDLAG ESEIHIPDVE PAAFLILLKY MYSDEIELEA DTVLATLYAA KKYIVPALAK ACVTFLETSL EAKNACVLLS QSRLFEEPEL TLRCWEVIDA QAELALHSEG FCEIDLQTL EILKRETLNT REAVVFQAAL DWAVAECKRQ GLGPTARNKR AVL GKALYLV RIPTMTLEEF ANGAAQSDVL TLEETHDVFL WYTAANKPKL EFPLQKRKGL TPQRCHRFQS SAYRSNQWRY RGRCDSIQFA VDKRIFIAGL GLYGSSGGKA EYSVKIELKR QGVTLAQNLT KFISDGSSNT FSVWFEHPVQ VEQDTFYTVS AVLDGNELSY FGQEGMTEVQ CGKVTFQFQC SSDSTNGTGV QGGQIPELVF YA
Source	Baculovirus
Target Names	btbd6b
Protein Names	Recommended name: BTB/POZ domain-containing protein 6-B
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.