



Recombinant Rat DDB1- and CUL4-associated factor 8 (Dcaf8)

Product Code	CSB-MP725920RA
Abbreviation	Dcaf8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U2M6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MSSKRPSTDG RRDLANGSLS SSPEEMSGAE EGRETSSGIE VEASDLSLSL TGDDGGPNRT STESRGTDTE SSGEEKDSDS MEDTGHYSIN DENRVHGHSD EEEEEEQPRH RVQRKRASRD QDSSDDERL EDWVSSETTA LPRPRWQALP ALRERELGSS ARFVYEACGA RVFVQRFRLQ HGLEGHTGCV NTLHFNQRGT WLASGSDDLK VVWDWVRRQ PVLDFESGHK SNVFQAKFLP NSGDSTLAMC ARDGQVRVAE LSATQCCKNT KRVAQHKGAS HKLALEPDSP CTFLSAGEDA VVFTIDLRQD RPASKLVVTK EKEKKVGLYT IYVNPANTHQ FAVGGRDQFV RIYDQRKIDE NENNGVLKKF CPHHLVNSES KANITCLVYS HDGTELLASY NDEDIYLFNS SHSDGAQYIK RYKGHRNNAT VKGVNFYGP K SEFVVS GSDC GHIFLWEKSS CQIIQFMEGD KGGVNCLEP HPHLPVLATS GLDHDVKIWA PTAEASTELT GLKDVIKKNK RERDEDSLHH TDLFDSHMLW FLMHHLRQRR HRRRWREPGV GATDADSDES PSSSDTSDEE EGPDRVQCMP S
Source	Mammalian cell
Target Names	Dcaf8
Protein Names	Recommended name: DDB1- and CUL4-associated factor 8 Alternative name(s): WD repeat-containing protein 42A
Expression Region	1-591
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.