



Recombinant Rat Pre-mRNA 3'-end-processing factor FIP1 (Fip111)

Product Code	CSB-MP725928RA
Abbreviation	Fip111
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U317
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MSAGEVERLV ELSSGGTGGDE EEEWLYGGPW DVHVHSDLAK DLDNEVERP EEENASANPP SGIEEEAAEN GVAKPKVTET EDDSDSDSDDD DEDDVHVTIG DIKTGAPQYG SYGTAPVNLN IKAGGRVYGN TGTKVKGVDL DAPGSINGVP LLEVLDLSFE DKPWRKPGAD LSDYFNNGFN EDTWKAYCEK QKRIRMGLEV IPVTSTTNKI TVQQGRTGNS EKEAALPSTK AEFTSPPSLF KTGLPPSRRL PGAIDVIGQT ITISRVEGRR RANENSNIQV LSDRSATEVD NNFSKPPPPF PPGAPPTHLP PPPFLPPPPT VSTAPPLIPP PGFPPPPGAP PPSLIPTIES GHSSGYDSRS ARAFPYGNVA FPHLTSSAPS WPSLVDTTKQ WDYYARREKD RDRDRERDRD RERERDRDRE RERTRERERE RDHSPTPSVF NSDEERYRYR EYAERGYERH RASREKEERH RERRHREKEE TRHKSSRSNS RRRHESEEGD SHRRHKHKKS KRSKEGKEAG SEPVPEQUEST EAAPAE
Source	Mammalian cell
Target Names	Fip111
Protein Names	Recommended name: Pre-mRNA 3'-end-processing factor FIP1 Alternative name(s): FIP1-like 1 protein
Expression Region	1-536
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.