



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R727 (MIMI_R727)

Product Code	CSB-EP725956ADAZ
Abbreviation	MIMI_R727
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UNY5
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MNNMDYDPRS MGSYGPNNYNN FNPNFGKIRR HKNTSYPGYN GFDPSDGPMN GPMGGPMGGP MNGPMNGQMG GPMNGSMNGP MNGPMNGPMN GQMGGPMNGP MGGPINGPMN GPRGRQMNGP NNGPMGGPMN GPMNGPNNNQ FNGPMNGPNE YYSPEDSDGS DYSDSNPNEF DSDDDDIDLS YFQKRQYDKP TYHKFEVTDE YLADYLLRKF GLDKDKVVRS IQKNKKDEFK AMILNYFYKE NPKIKHQSIN NQYKLFKRKWN KSGLKIKIDD LETYYENKVK PHMH
Source	E.coli
Target Names	MIMI_R727
Protein Names	Recommended name: Uncharacterized protein R727
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.