



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R588 (MIMI_R588)

Product Code	CSB-YP725961ADAZ
Abbreviation	MIMI_R588
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UP50
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MSSFIFPTYR EVHSLQSTVS QLASNSGVKN KLANYGLNVS SVAWEDTSRY LNSCVGSNIS DVTLKVNKSR MPIIRSNNFT DVTVDVPVSK LPKLVVGNHD GSQKTVTLE EYLKDFHLYC GEFTRDRYNL FADRDQHVLT SSQACVLPLN SGEVEFAVDL YNYQSRNEPA VLVIVSTAYG TSAQVVCGGN TELYFNDNSV SRMFKAERIE DFRKSQGRST SGPMTTEEKA LNGIYIYQVP LVIKQERPNI ESCFLGFVEC DSLEASATFS NDSNNFFEHS ARNSRGMSRA ILSLGQEKGP YKGIKKTGG IYKLERDQTK PIRLTVQFYM CTDITNLSDN NVEEISNQIR RIYDQGLNEG SLVVDVSTTK PGLHQPTQKR PTQTTSKPVI NDEITML
Source	Yeast
Target Names	MIMI_R588
Protein Names	Recommended name: Uncharacterized protein R588
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.