



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R839 (MIMI_R839)

Product Code	CSB-YP726022ADAZ
Abbreviation	MIMI_R839
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQI6
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MELLEKFISQ DINNNFLHYL NFKDNGSKKY ITNQLIFLEK LFDIPYGNLF LFTKMF DNKV WEIPTNELSE GLVRLFDRLE LTKINELAAG NGLLSARLKY YSEKLNTNLK IKTS DGSSKN FGNHEFTYIK VSESDISDFD KSGPIISWI HRFFEHELLS YVKKHNN EYI FLIGEHPDST DYGN NHSMYF HNKMY SFGYQ HQIIEFQQIS QMDYYCYDKI RNDIYNENKT CVVLYYKTHL HSKVTNIINT LKQNH PKLFG KFLNKNKKYY DQDKILLEIS KNKINEYFSY CHDDLCLLLT DRYKKLVES LMDNDDSDSE ISHKTIQLSN IIKMMYCGLT KVNLRPIIS CSISNCNSEQ SLSRIVRHNE NKPTTIILNS VKNNNSVISH SVDDKLSQRY GKKGTIGFSI KRQEYEQENK YTSPFINGSS IPTRMIIRDI TERIFAEIMH IKSNINI
Source	Yeast
Target Names	MIMI_R839
Protein Names	Recommended name: Uncharacterized protein R839
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.