



Recombinant Acanthamoeba polyphaga mimivirus Putative FNIP repeat-containing protein L415 (MIMI_L415)

Product Code	CSB-MP726025ADAZ
Abbreviation	MIMI_L415
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQL6
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MNIVEMLCDT VIVYMMDYLS DIDKINFLKT CSSFYNHRFC VKFNDFHSFD KIRDLDYSEN FKRVFTKTK NNVKSFTLNN SYDRIPIDAK HLIIGKLFKL SLENIIPRGI THLTIEEEYG SEINPGIINP EEIPDIVTYL NIKSRNIFIK KGAIPDSVTH LYFGSDYLSK DIIPKNVVYL RFGDFSRCVI GQRLVETRQV NDITVTTTRI ESYIPESVAC LHLGFMSMSNK CLGERIETQ EMDGEKINTT KIVSYIPKNV TCLTIFNEVY KNTDFDFIPKN VTKLNLFG EI MEMIDETIIP KNITTLGLAC YRDLNKIKIP KKVTRLNFLF EEYMSNPVKL DIPEHITTIS LKYTHDIFKK VVFNNILNLK IDSYYFEKKL YLPNSIKTLT IYGN YDKKVI KKLPNTINTL IFKHK
Source	Mammalian cell
Target Names	MIMI_L415
Protein Names	Recommended name: Putative FNIP repeat-containing protein L415
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.