



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R326 (MIMI_R326)

Product Code	CSB-BP726033ADAZ
Abbreviation	MIMI_R326
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQR5
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MEDPIKIIHK YKNNNGRIQY HINIFIGDIV DENCMRILRK IKNLDLYTSL TSLETREIDI LEKNYGEYWY EKFFNSYHIN NTKELTLKNS VRMRELRSLY TEEWVNRHFV NYKKRLETTI FNYEYVVKED RERRSVKRHI RRQQDDVEEL LDYRTTGHP LPPALNDEQYV MSRIESSVET DNWCADNLLE EQLTESELNR ELKKLVSKQN LDYSTDKPED SESEDIELED SESESESED IDQHGGQGPD DDEFNANFDD PQFDEFDFGN TEDEQNVKLF EMEVEQDLED VDLLFNDIDE TDKNSKLTTR EIKEAISNEQ YDRIGKKIVD FDES RDNSMF DENLKDVITK TYITNQYLFK DDTIKTIRDK ICCGFKNSNK FGENTYIIPS HQYLVSEYTY QGKIDRVMIG QKWIIRNDIL KLDVEPNTNT SVKKNSGEIY VF
Source	Baculovirus
Target Names	MIMI_R326
Protein Names	Recommended name: Uncharacterized protein R326
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.