



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R369 (MIMI_R369)

Product Code	CSB-MP726038ADAZ
Abbreviation	MIMI_R369
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQV5
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MTDNINIVHK KLFNVIDANF TKILKLKGVV NDIDDNGKFL TIKTLDTERT LKCRTKDNFT KVDVGDVVDV VGNLQLDSND LSNIYLDLKY LSKQSTFDPS KALSIHNKLL VYLNEHEKVS KVKKIHKKS LPKQIYNVAL IAFPTGINLV NTFKTEFKNK CVGNLYTFYM DTDKTDVSIV SALEYFRKYR NIDIICLLTD QLNIMETCKL SSKEIAKVFL SRKNCPYILS IVSDSSEKNP EPITALLSNK KIETTLAAVN FISEVQWEYR QKVNRAVDLV KEKFYIHMDK KIKKLQQYEV LIWNLNKTSP PINSTDNLKN LLRNKMTEEF DRLNRIETEL VNNIIMSDPI QKIVNNIDMF DSNIFMKIVD SIIDKMKNDK AQDKSLATTI KENALRQNLA FQSGQVDINK SRTLQDALVE EMKD HITQDP YLNYKERDEN DIVDP
Source	Mammalian cell
Target Names	MIMI_R369
Protein Names	Recommended name: Uncharacterized protein R369
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.