



# Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R369 (MIMI\_R369)

<b>Product Code</b>	CSB-YP726038ADAZ
<b>Abbreviation</b>	MIMI_R369
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5UQV5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acanthamoeba polyphaga mimivirus (APMV)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTDNINIVHK KLFNVIDANF TKILKLGVV NDIDDNGKFL TIKTLDTERT LKCRTKDNFT KVDVGDVVDV VGNLQLDSND LSNIYLDLKY LSKQSTFDPS KALSIHNKLL VYLNEHEKVS KVVKKIHKKS LPKQIYNVAL IAFPTGINLV NTFKTEFKNK CVGNLYTFYM DTDKTDVSIV SALEYFRKYR NIDIICLLTD QLNIMETCKL SSKEIAKVFL SRKNCPYILS IVSDSSEKNP EPITALLSNK KIETTLAAVN FISEVQWEYR QKVNRAVDLV KEKFYIHMDK KIKKLQQYEV LIWNLNKTSP PINSTDNLKN LLRNKMTEEF DRLNRIETEL VNNIIMSDPI QKIVNNIDMF DSNIFMKIVD SIIDKMKNDK AQDKSLATTI KENALRQNLA FQSGQVDINK SRTLQDALVE EMKD HITQDP YLNYKERDEN DIVDP
<b>Source</b>	Yeast
<b>Target Names</b>	MIMI_R369
<b>Protein Names</b>	Recommended name: Uncharacterized protein R369
<b>Expression Region</b>	1-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.