



# Recombinant Acanthamoeba polyphaga mimivirus Putative ankyrin repeat protein R580 (MIMI\_R580)

<b>Product Code</b>	CSB-BP726051ADAZ
<b>Abbreviation</b>	MIMI_R580
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5UR56
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acanthamoeba polyphaga mimivirus (APMV)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSKTYIIST NDYFDPINLS IKDGGEKTID IIDCERYTLI DCIIFFCKKE SFVIEIIVPD FLTCKDLLCD SENLSLCVVN NIYYMNDISV FDSMSLEKST LKSLMNWAVK NNYKVLLQKY LDNQNIKIYD LLACGISNSN VCDNILELLL SYCKYFQENE TMFLVSFCLE MHRIDIATL IMYNSYPSVI NKSLNFASKS NMSELAIFLV DNGAEINCDH ILCFIRIGDI DTVLYMLCRY DVQKLIDRCH FDLVEAVFNS GSLEMIETFI DFGMKINSKV YQHLKYDSKD TYEILRILLA REIYPRKCSV VLKNAAHICS IDVLKILVDF EFNQKSFDNA LVSTINAGKF KNAEYLLFSG ANINFNMPPT NCMFKINFQT IKFLIDNNFD LEIHGTLILN KSLNNGYYDC ANILIENGVK FSLTKSELLK IYSDFQNTTH GQIDPELNIN SLTSDELENL IYIDIVNK
<b>Source</b>	Baculovirus
<b>Target Names</b>	MIMI_R580
<b>Protein Names</b>	Recommended name: Putative ankyrin repeat protein R580
<b>Expression Region</b>	1-468
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.