



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized hydrolase R843 (MIMI_R843)

Product Code	CSB-YP726060ADAZ
Abbreviation	MIMI_R843
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5URB6
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MDNDDYIWLD YLTVEQPNQE KTDISVKHYY GIDIKTDKYL QNAIDIWELI NKYGKKYFCN YGGETFQINI GDRFYVAHFP SQCLNNGCKY PVLIFLHGLT GYSWHSALDK TGLIELANQN NFIVLFGQGN GEITRPVRDK YGGVSFGDIY WQIENPELDM DYIRNVMLNQG GTIEYCRQDD LNLLELRKMF SDKIYLGYS NGAMYSFNMCFYDLHFSGIC SMMGGYGGGLA GYSNSKISEI VDRPINISLD IPIILSGSL DEYLPASKKA FDILLGKNFS NVKFLSIPDR KHTYSRDFEK YIWEFFFTST KN
Source	Yeast
Target Names	MIMI_R843
Protein Names	Recommended name: Uncharacterized hydrolase R843 EC= 3.-.-
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.