



Recombinant *Drosophila simulans* Pro-corazonin (Crz)

Product Code	CSB-BP726130DMJ
Abbreviation	Crz
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5W1L5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila simulans</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	Q TFQYSRGWTN GKRSFNAASP LLANGHLHRG SELGLTDLYD LQDWSSDRRL ERCLSQLQRS LIARNCVPGS DFNANRVDPD PENSVHPRLS NINGENVLYS SANIPNRHRQ SNELLEELSA AGGASAEPNV FGKH
Source	Baculovirus
Target Names	Crz
Protein Names	Recommended name: Pro-corazonin Short name= Crz Short name= Ds-Crz Cleaved into the following 2 chains: 1. Corazonin 2. Corazonin precursor-related peptide Short name= 3. CPRP
Expression Region	20-154
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.