



# Recombinant *Bacillus clausii* Ribose import ATP-binding protein RbsA (rbsA)

<b>Product Code</b>	CSB-MP726140BAAC
<b>Abbreviation</b>	rbsA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5WC31
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Bacillus clausii</i> (strain KSM-K16)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MIVEMSQISK SFSGNQVLKN VSFTLEKGEI HALMGENGAG KSTLMKILTG IYPRDEGVVR VKGKEVAYTH PKDAEKDGIA VIHQELNILP ELSVAENLFL GNEYTVGKTG WLKTKKMNEY AEQQLAELGL KVKASDRAGA LSVGKQQLIE IAKALMTNAD IIIMDEPTAA LTDREIDTLF TVIRGLQEKG VTFVYISHRM EEIFSLCQRI TVLRDGNYS TKVIAETSD EIVKLMVGRE LGGRYPDHDL EPGAVKLAVK NVSRKNEFES VTFEVREGEI FGIAGLMGAG RSEVVESIFG YRHLHSGTIE IDQKPVVKR PVDIENGIA FVPEDRKTG LIVNFSVRDN LSITNLATVS SGNMMSSTKE QALYQELVKK LGVRTSGPHQ KAKSLSGGNQ QKVVIKWLK MKPKVLILDE PTRGVDVGAK KEIYTIMNEL AKQGVAIMV SSELPEIIGM STRVAVMFEG KLMKVLDRNE LSEETIMHYA TGGDKHVHQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rbsA
<b>Protein Names</b>	Recommended name: Ribose import ATP-binding protein RbsA EC= 3.6.3.17
<b>Expression Region</b>	1-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.