



Recombinant *Bacillus clausii* Hydroxyethylthiazole kinase (thiM)

Product Code	CSB-MP726146BAAC
Abbreviation	thiM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5WDW4
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus clausii</i> (strain KSM-K16)
Purity	≥85% (SDS-PAGE)
Sequence	MELVRQNNPL IHCMTNDVVM NFTANGLLAI GASPVMASFSA AETEDMASVA DALLLNIGTP SNEGVDSMVL AGKAANRAGK PVIFDPVGVG ATPFRNEISK RILNEVDVAV VRGNAGEIAA LLGQAGTVKG VDGKINADVK ELCLQAAKAL RTVVVVTGEV DAVSNGEQLV AVANGHEWLT KVVGTGCLLG GVIAAYA AVQ PNSLVDAAVE ALAFYGVCAE QAYEQTKKEG IGSFQQTFLN ELGNMTDEKA AALKRVEQIC
Source	Mammalian cell
Target Names	thiM
Protein Names	Recommended name: Hydroxyethylthiazole kinase EC= 2.7.1.50 Alternative name(s): 4-methyl-5-beta-hydroxyethylthiazole kinase Short name= TH kinase Short name= Thz kinase
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.