



Recombinant Legionella pneumophila Elongation factor P (efp)

Product Code	CSB-YP726242LAAA
Abbreviation	efp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5WZP1
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila (strain Lens)
Purity	≥85% (SDS-PAGE)
Sequence	MAVYSTNEFK NGLKVMVDDA PCSILDCEFV KPGKGQAFTR IKIRNLKTGR VVERTFKSGD TLPSADVADV EMQYLYNDGE HWHFMVPDTF EQYAVTENIL ADAAQWLKEQ DVCVVTWNN EPIQVTPPNF VILAITETDP GLKGDTSGGG GKPATLETGA VVRVPLFVQT GELIKVDTRK GEYVSRAKE
Source	Yeast
Target Names	efp
Protein Names	Recommended name: Elongation factor P Short name= EF-P
Expression Region	1-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.