



Recombinant Tryptophan--tRNA ligase (trpS)

Product Code	CSB-EP726290SMZ-B
Abbreviation	trpS
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5X9A2
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MTKPIILTGD RPTGKLHLGH YVGLKNRVF LQENENKYKMF VFLADQQALT DHAKESLIQ ESIGNVALDY LSVGLDPKQS TIFIQSQIPE LAELSMYYMN LVSLARLERN PTVKTEIAQK GFGESIPSGF LVYVPVQAAD ITAFKANLVP VGNDQKPMIE QTREIVRSFN HTYHTDCLVE PEGIYPENEK AGRLPGLDGN AKMSKSLGNG IYLSDDADTV RKKVMSMYTD PNHKIEDPG QIEGNMVFHY LDIFARKEDQ ADIEAMKEHY QRGGLGDVKT KRYLLDILER ELAPIRERRL EYAKDMGEVF RMLQEGSQKA RTVAAKTLSE VKSAMGINYF
Source	E.coli
Target Names	trpS
Protein Names	Recommended name: Tryptophan--tRNA ligase EC= 6.1.1.2 Alternative name(s): Tryptophanyl-tRNA synthetase Short name= TrpRS
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>