



# Recombinant UDP-glucose 6-dehydrogenase (hasB)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP726291SMZ-B                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Abbreviation</b>      | hasB                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                              |
| <b>Uniprot No.</b>       | Q5X9A8                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Sequence</b>          | MKIAVAGSGY VGLSLGVLLS LQNEVTIVDI LPSKVDKINN GLSPIQDEYI<br>EYYLKSKQLS IKATLDSKAA YKEAELVIA TPTNYNSRIN YFDTQHVETV<br>IKEVLSVNSH ATLIKSTIP IGFITEMRQK FQTDRIIFSP EFLRESKALY<br>DNLYPRIIV SCEENDSPKV KADAEKFALL LKSAACKNNV PVLIMGASEA<br>EAVKLFANTY LALRVAYFNE LDTYAESRKL NSHMIIQGIS YDDRIGMHYN<br>NPSFGYGGYC LPKDTKQLLA NYNNIPQTLI EAIVSSNNVR KSYIAKQIIN<br>VLKEQESPVK VVG VYRLIMK SNSDNFRESA IKDVIDILKS KDIKIIYEP<br>MLNKLESEDQ SVLVNDLENF KKQANIIVTN RYDNELQDVK NKVYSRDIFG<br>RD |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Target Names</b>      | hasB                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | Recommended name: UDP-glucose 6-dehydrogenase Short name= UDP-Glc dehydrogenase Short name= UDP-GlcDH Short name= UDPGDH EC= 1.1.1.22                                                                                                                                                                                                                                                                                                                                           |
| <b>Expression Region</b> | 1-402                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                     |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                               |



of lyophilized form is 12 months at -20°C/-80°C.