



# Recombinant Probable dipeptidase B (pepDB)

<b>Product Code</b>	CSB-EP726297SMZ-B
<b>Abbreviation</b>	pepDB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5X9M0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MINKKISLGV LSILTAFLSQ SVSYACTGFI IGKDLTKDGS LLXGRTEDE PHHNKNFIVR LAKDNPAGEK WKDLSNGFEY PLPEHSYRYS AIPDVTPNKG VYDEAGFNEC GVSMSATVSA SANDAIQKID PYVKNGLAES SMASVILPSV KTAREGVALI AKIVTEKGAA EGNIVTLADK DGIWYMEILS GHQYVAIKFP DDKYAVFPNT FYLGHVDFND KENTIASADV EKVAKKAKSY IEVDGKFHIA KSYNPPLNDA NRSRSFSGIK SLDPDSKVTY KDSNYELLQS TDKTFSLEDA MKLQRNRFEG LDLKPLDQMA LDGKGKPKSK KAVKGYAYPI SNPNVMEAH FQLKKDIPAE LGGGVMWLSI GSPRNAPYLP YLGNISRTYE AYQEKSTQYN DKSWYWTVSH INDLVAAHPK PFGTKVIDEI KGLEKTWIAE QDKTTKEISD LVVSDPKAAQ EKADKISLDR AEKTFKRLKA IEAKLVKEKP KNKKGLNRS
<b>Source</b>	E.coli
<b>Target Names</b>	pepDB
<b>Protein Names</b>	Recommended name: Probable dipeptidase B EC= 3.4.-.-
<b>Expression Region</b>	1-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.