



# Recombinant Asparagine--tRNA ligase (asnS)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP726332SMZ   |
| <b>Abbreviation</b>      | asnS  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q5XD20  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSKKLISIVN VKDYVGQEV T IGAWVANKSG KGKIAFVQLR DGSAFFQGVA<br>FKPNFIEKYG EESGLEKFDV IKRLNQETSV YVTGIVKEDE RSKFGYELDI<br>TDLEIIGESH EYPVTPKEHG TDFLMDNRHL WLSRKRQMAV MQRNIIYA<br>TYEFFDQNGF IKFDSPILSE NAAEDSTELF ETDYFGKPAF LSQSGQLYLE<br>AGAMALGRVF DFGPVFRAEK SKTRRHLTEF WMMDAEYSFL SHEESLDLQE<br>AYVKALIQQV LDRAPQALDI LERDVEALKR YITEPFKRVS YDDAITLLQE<br>HEADEDTDYE HLEHGDDFGS PHETWISNYF GVPTFVVNYP ASFKAFYMKP<br>VPGNPERVLC ADLLAPEGYG EIIGGSMRED NYDALVAKMD ELGMDKSEYD<br>FYLDLRKYGS VPHGGFGIGI ERMVTFVAGT KHIREAIPFP RMLHRIRP |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | asnS  |
| <b>Protein Names</b>     | Recommended name: Asparagine--tRNA ligase EC= 6.1.1.22 Alternative name(s): Asparaginyl-tRNA synthetase Short name= AsnRS   |
| <b>Expression Region</b> | 1-448   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |