



Recombinant Asparagine--tRNA ligase (asnS)

Product Code	CSB-MP726332SMZ
Abbreviation	asnS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XD20
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MSKKLISIVN VKDYVGQEV T IGAWVANKSG KGKIAFVQLR DGSAFFQGVA FKPNFIEKYG EESGLEKFDV IKRLNQETSV YVTGIVKEDE RSKFGYELDI TDLEIIGESH EYPVTPKEHG TDFLMDNRHL WLSRKRQMAV MQRNIIYA TYEFFDQNGF IKFDSPILSE NAAEDSTELF ETDYFGKPAF LSQSGQLYLE AGAMALGRVF DFGPVFRAEK SKTRRHLTEF WMMDAEYSFL SHEESLDLQE AYVKALIQGV LDRAPQALDI LERDVEALKR YITEPFKRVS YDDAITLLQE HEADED TDYE HLEHGDDFGS PHETWISNYF GVPTFVVNYP ASFKAFYMKP VPGNPERVLC ADLLAPEGYG EIIGGSMRED NYDALVAKMD ELGMDKSEYD FYLDLRKYGS VPHGGFGIGI ERMVTFVAGT KHIREAIPFP RMLHRIRP
Source	Mammalian cell
Target Names	asnS
Protein Names	Recommended name: Asparagine--tRNA ligase EC= 6.1.1.22 Alternative name(s): Asparaginyl-tRNA synthetase Short name= AsnRS
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.