



Recombinant Arabidopsis thaliana Methyl-CpG-binding domain-containing protein 1 (MBD1)

Product Code	CSB-MP726350DOA
Abbreviation	MBD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XEN5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MLPFPAMNLK KSRSENSVA SSGSKIEEQT EKSAEPTTIK VQKKAGTPGR SIDVFAVQCE KCMKWRKIDT QDEYEDIRSR VQEDPFFCKT KEGVSCEDVG DLNYDSSRTW VIDKPLPRT PRGFKRSLIL RKDYSKMDAY YITPTGKKLK SRNEIAAFID ANQDYKYALL GDFNFTVPKV MEETVPSGIL SDRTPKPSRK VTID
Source	Mammalian cell
Target Names	MBD1
Protein Names	Recommended name: Methyl-CpG-binding domain-containing protein 1 Short name= AtMBD1 Short name= MBD01 Alternative name(s): Methyl-CpG-binding protein MBD1
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.