



# Recombinant *Xenopus tropicalis* Ribose-phosphate pyrophosphokinase 2 (prps2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP726356XBF                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Abbreviation</b>      | prps2                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q5XGI0                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | <i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | PNIVLFSGS SHQDLSQKVA ERLGLELGKV VTKKFSNQET SVEIGESVRG<br>EDVYIIQSGC GEINDNLMEL LIMINACKIA SSSRVTAVIP CFPYARQDKK<br>DKSRAPISAK LVANMLSVAG ADHIITMDLH ASQIQGFFDI PVDNLYAEPA<br>VLQWIRENIP EWKNSIIVSP DAGGAKRVTS IADRLNVEFA LIHKERKKAN<br>EVD RMVLVGD VKDRVAILVD DMADTCGTIC HAADKLLSAG ATKVYAILTH<br>GIFSGPAISR INNAAFEAVV VTNTIPQEDK MKQCCKSIQVI DISMILAEAI<br>RRTHNGESVS YLF SHVPL                           |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | prps2                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Ribose-phosphate pyrophosphokinase 2 EC= 2.7.6.1<br>Alternative name(s): Phosphoribosyl pyrophosphate synthase II Short name= PRS-II                                                                                                                                                                                                                                                      |
| <b>Expression Region</b> | 2-318                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |