



Recombinant Danio rerio Protein TSSC1 (tssc1)

Product Code	CSB-YP726387DIL
Abbreviation	tssc1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5XJP1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEDDAPVIYG LEFQARALTA QTAETDAIRF LVGTQSLKFD NQIHIDFDD ENNMINKSVL LHDAGEIWHI SGSPADKAVL TTCYNKTSES RVLTCAAVWR MPPEWESGNH ESPDDSSNNP QTLELLCHLD NSAHGSMACV LWPEMGDGKR VISLAENHVL LWDLQESSAK ATISSSATLE GKGQLKFTAG KWSPHHNCTQ LATANDTAIR GWDLRMSQI YCIENAHGQL VRDLDFNPKNK QYYLASCDD CKVKFWDVRH ISEPVKCLEE HSHWVWSVRY NHSHDQLLLT ASSDSRVILS NMVSISSSEPF GHLVDDDDLS DPEENQQEDK GKEPLQDSII STYEEHEDSV YAVEWSAADP WLFASLSYDG RLVINRV PRA LKYRILL
Source	Yeast
Target Names	eipr1
Protein Names	Recommended name: Protein TSSC1
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>