



Recombinant Ceratonia siliqua Maturase K (matK)

Product Code	CSB-MP726410CGAT
Abbreviation	matK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5YK50
Product Type	Recombinant Protein
Immunogen Species	Ceratonia siliqua (Carob) (St John's-bread)
Purity	>85% (SDS-PAGE)
Sequence	MEEFQVYLEL ARSRQHDFLY PLIFREYIYA LAYDHGLNSS ILVENLGYDN KSSLLIVKRL ITRMYQQNHL IISANDSNKK RFWGYNKNLY SQIISEGFAV SVEIPFSLQL ISSLEEAKII KSYNLRISIHS IFPFEDKFP YLNYVSDVRI PYPIHLEILV QILRYWVKDA SSFHLLRLFY YEYYNWNLSLI TPKKWISTFS KSNPRFFLFL YNFYVCEYES IFLFLRNKSS YLRLTSSGVF FERIYFYAKI EHLVEVFDKD FPSTLWFFKD PFIHYVRYQG KSILASKNTP FFMNKWKYYL IHLWQCNFYV WSQPGKIHIN QLSEHSFYFL GYFSNURLNP SVVRSQMLEN SFIENVMKK LDTIPIIPL IRSLAKVKFC NGLGHPISKP VWADSSDFDI IDRFLRICRN LSHYYNGSSK KKSlyRIKYI LRLSCIKTLA RKHKSTVRVF LKRLGSELLE ELFTEEEEIL SLIFPRASST LHRLYRGRIW YLDIFYFHQ
Source	Mammalian cell
Target Names	matK
Protein Names	Recommended name: Maturase K Alternative name(s): Intron maturase
Expression Region	1-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.