



Recombinant Chicken Serine/threonine-protein phosphatase 2A 55 kDa regulatory subunit B delta isoform (PPP2R2D)

Product Code	CSB-EP726492CH
Abbreviation	PPP2R2D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZIY5
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MAGVAGGGGN GNDFQWCFSQ VQGAVDEDVA EADIISTVEF NYSGDLLATG DKGGRVVIFQ REQENKSRPH SRGEYNVYST FQSHEPEFDY LKSLEIEEKI NKIRWLPQQN AAHFLSTND KTIKLWKISE RDKRAEGYNL KDEDGRLRDP FKITALRVPI LKPMDLMVEA SPRRIFANAH TYHINSISVN SDHETYLSAD DLRINLWHLE ITDRSFNIVD IKPANMEELT EVITAAEFHP HHCNVFVYSS SKGTIRLCDM RSSALCDRHS KFFEEPEDPS SRSFFSEIIS SISDVKFSHS GRYMMTRDYL SVKVWDLNME NRPVETYQVH EYLRSKLCSL YENDCIFDKF ECCWNGSDGA IMTGSYNFF RMFDRNTRRD ITLEASRESS KPRAILKPRK VCTGGKRKKD EINVDSLDFN KKILHTAWHP MENIIA VAAT NNLYIFQDKI N
Source	E.coli
Target Names	PPP2R2D
Protein Names	Recommended name: Serine/threonine-protein phosphatase 2A 55 kDa regulatory subunit B delta isoform Alternative name(s): PP2A subunit B isoform B55-delta PP2A subunit B isoform PR55-delta PP2A subunit B isoform R2-delta PP2A subunit B
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.