



Recombinant Chicken Anamorsin (CIAPIN1)

Product Code	CSB-MP726497CH
Abbreviation	CIAPIN1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5ZJA2
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MGEYGIAPGQ RVAVIWDSSS PVEALKGLVD AVQASVGADS RVSVENINQL CQSAHRESSF DVILSGVVPD STAQHSAEVL AEIARILKPG GRVLLKEPVV TESENNSIK TAAKLPALT LSGLVEVKGL QKEPLTAEAA QSVREHLGYQ GNDLLIVQIE GRKPNFEVGS SSQKLSFAK KTSPSGKPSV DPATAKLWTL SASDMNDEEM DLLDSELDL SEDLKKPDPA SLRAPSCKEK GKKKACKNCT CGLAEELEQE KKSSQPKSAC GNCYLGDAFR CASCPYLGDA CLQAWREDPA EREPAA
Source	Mammalian cell
Target Names	CIAPIN1
Protein Names	Recommended name: Anamorsin Alternative name(s): Cytokine-induced apoptosis inhibitor 1 Fe-S cluster assembly protein DRE2 homolog
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>