



Recombinant Chicken Cytosolic Fe-S cluster assembly factor NUBP2 (NUBP2)

Product Code	CSB-BP726524CH
Abbreviation	NUBP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZKV4
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MEEVVAERSN LGGVRHILLV LSGKGGVGKS TISTELALSL RHSGKKVGIL DVDLCGPSIP RMFKVQDNDV HQCDAGWVPV FVDQEKSISL MSIGFLLEKP DDAVVWRGPK KNALIKQFVA DVAWGELDFL IVDTPPGTSD EHISTVEALR PYKPLGAILV TTPQAVSVGD VRRELTFCCK TGLRVLGIVE NMSGFVCPHC SECTNIFSKG GGEELAKHAG VPFLGSVPLD PQLSQSLEEG RDFIQEFPKS SAFPALTRIA QQILDGALQR SS
Source	Baculovirus
Target Names	NUBP2
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NUBP2 Alternative name(s): Nucleotide-binding protein 2 Short name= NBP 2
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.