



Recombinant Chicken Protein FAM46C (FAM46C)

Product Code	CSB-EP726529CH
Abbreviation	FAM46C
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5ZL95
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	<p>MASKGSNTDC MSCSVLNWEQ VSRLHEVLTE VVPIHGRGNF PTLKITLKDI VQTVRSRLNE AGIAVHDVRL NGSAAAGHVLV KDNGLGCKDL DLIFQVSLPS EAEFQLVRDV VLRSLNFLP EGVSKLKISP VTLKEAYIQK LVKVSTESDR WSLISLSNKH GRNVELKFVD CIRRQFEFSV DSFQIILDSL LFYYDYSENP MSEHFHPTVI GESMYGDFA AFDHLQNKLI ATKNPEEIRG GLLKYSNLL VRDFRPMDKD EIKTLERYMC SRFFIDFPDI LDQQRKLETY LQNHFSSKEER SKYDYLMILR RVVNESTVCL MGHERRQTLN LISLLALKVL AEQNIIPNAT TVTCYYQPAP YVSDVNFSNY YLANPPVPYS QSYPTWLPCN</p>
Source	E.coli
Target Names	FAM46C
Protein Names	Recommended name: Protein FAM46C
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>