



Recombinant Chicken ATP synthase subunit beta, mitochondrial (ATP5B)

Product Code	CSB-EP726530CH-B
Abbreviation	ATP5B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZLC5
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	QAAPAAK AGSATGRIVA VIGAVVDVQF DEGLPPILNA LEVQGRETRL VLEVAQHLGE NTVRTIAMDG TEGLVRGQKV LDSGAPIRIP VGPETLGRIM NVIGEPIDER GPITTKQFAA IHAEAPEFVE MSVEQKILVT GIKVVDLLAP YAKGGKIGLF GGAGVGKTVL IMELINNVAK AHGGYSVFAG VGERTREGND LYHEMIESGV INLKDATSKV ALVYGQMNEP PGARARVALT GLTVAEYFRD QEGQDVLLFI DNIFRFTQAG SEVSALLGRI PSAVGYQPTL ATDMGTMQER ITTRKGSIT SVQAIYVPAD DLTPAPATT FAHLDATTVL SRAIAELGIY PAVDPLDSTS RIMDPNIVGP EHYDVARGVQ KILQDYKSLQ DIIAILGMDE LSEEDKLTVA RARKIQRFLS QPFQVAEVFT GHMGKLVPLK ETIKGFKQIL AGEYDHLPEQ AFYMGVPIEE AVAKAEKLA EHA
Source	E.coli
Target Names	ATP5B
Protein Names	Recommended name: ATP synthase subunit beta, mitochondrial EC= 3.6.3.14
Expression Region	54-533
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a subunit of mitochondrial ATP synthase. Mitochondrial ATP synthase catalyzes ATP synthesis, utilizing an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. ATP synthase is composed of two linked multi-subunit complexes: the soluble catalytic core, F1, and the membrane-spanning component, Fo, comprising the proton channel. The catalytic portion of mitochondrial ATP synthase consists of 5 different subunits (alpha, beta, gamma, delta, and epsilon) assembled with a stoichiometry of 3 alpha, 3 beta, and a single representative of the other 3. The proton channel consists of three main subunits (a, b, c). This gene encodes the



beta subunit of the catalytic core.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.