



Recombinant Chicken Transcription initiation factor TFIID subunit 8 (TAF8)

Product Code	CSB-YP726544CH
Abbreviation	TAF8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZMS1
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MLQSYISEIG RSAKSYCEHT ARTQPTLSDI VVTLVEMGFN VETLPAYAKR SQRMVITAPP VTNPQVTPKA LTAGQNKPHP SHIPGYFPEF PDPHTYIKTP TYREPVS DYQ VLREKAASQR RDVERALTRF MAKTGETQSL FKDDVSTFPL IAARPFTVPY LTALLPSELE MQQMEETDSS EQDDQTD TEN LPLHISTDDS GPEKENASVL QQNTSLSGSR NGEENMIDNP YLRPVKKPKI RRKK
Source	Yeast
Target Names	TAF8
Protein Names	Recommended name: Transcription initiation factor TFIID subunit 8 Alternative name(s): TBP-associated factor 8
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes one of several TATA-binding protein (TBP)-associated factors (TAFs), which are integral subunits of the general transcription factor complex TFIID. TFIID recognizes the core promoter of many genes and nucleates the assembly of a transcription preinitiation complex containing RNA polymerase II and other initiation factors. This protein contains an H4-like histone fold domain, and interacts with several subunits of TFIID including TBP and the histone-fold protein TAF10. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.