



Recombinant Legionella pneumophila subsp. pneumophila Phosphoglycerate kinase (pgk)

Product Code	CSB-EP726591LDL-B
Abbreviation	pgk
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZZ74
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)
Purity	≥85% (SDS-PAGE)
Sequence	MNLIKMSDID LSGKRVLIRE DLNVPIKDGM ITSDQRLQAA LPTIKSALDS GAAVIVLSHL GRPEEGKYEK KFSLEPVADY LRENLEYPVR FVKDYLSGVG VNPGEVLVCE NVRFNPGEKA NDEALAKKLA SLCDVFVMDA FGTAHRAQAS TYGVAQYAPV AVAGPLLIRE LEALNQVLKA PKKPIVAIVG GAKVSSKLSL LKQLVGMVDV LIPGGGIANT FLKAQGFEIG ISLYEPDLLD EARHILILAK EKGCQIPLPT DVVVGKTFSE TCPAFNKSLN NVAADDMLD IGPETIRDYV DLIHDANTII WNGPVGVFEE PQFAYGTRAI AIAIAESDAF SIAGGGDTLA AVDLYDLNQQ ISYISTGGGA FLECLEGKTL PAVAILQERA KHVKTN
Source	E.coli
Target Names	pgk
Protein Names	Recommended name: Phosphoglycerate kinase EC= 2.7.2.3
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.