



# Recombinant Mouse BAG family molecular chaperone regulator 1 (Bag1)

<b>Product Code</b>	CSB-MP726670MO
<b>Abbreviation</b>	Bag1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q60739
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGRSAARRP RGDREPLGPR LRAPRPAREP RQSESRAERG LPPSQRSSVR SAASGHDRST RGAPAGACKP RVKKKVRPRS SQSEKVGSSS RELTRSKKVT RSKNVTGTQV EEVTKIEEAT QTEEVTVAAE VTQTDNMAKT EEMVQTEEME TPRLSVIVTH SNERYDLLVT PQQGNSEPVV QDLAQLVEEA TGVPLPFQKL IFKGKSLKEM ETPLSALGMQ NGCRVMLIGE KSNPEEEVEL KKLKDLEVSA EKIANHLQEL NKELSGIQQG FLAKELQAEA LCKLDRKVK TIEQFMKILE EIDTMVLPEQ FKDSRLKRKN LVKKVQVFLA ECDTVEQYIC QETERLQSTN LALAE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Bag1
<b>Protein Names</b>	Recommended name: BAG family molecular chaperone regulator 1 Short name= BAG-1 Alternative name(s): Bcl-2-associated athanogene 1
<b>Expression Region</b>	1-355
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The oncogene BCL2 is a membrane protein that blocks a step in a pathway leading to apoptosis or programmed cell death. This protein binds to BCL2 and is referred to as BCL2-associated athanogene. It enhances the anti-apoptotic effects of BCL2 and represents a link between growth factor receptors and anti-apoptotic mechanisms. At least three protein isoforms are encoded by this mRNA through the use of a non-AUG (CUG) start site, and alternative, downstream, AUG translation initiation sites.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.