



# Recombinant Mouse Nuclear factor erythroid 2-related factor 2 (Nfe2l2)

<b>Product Code</b>	CSB-BP726674MO
<b>Abbreviation</b>	Nfe2l2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q60795
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMDLELPPPG LQSQQDMDLI DILWRQDIDL GVSREVFDFS QRQKDYELEK QKKLEKERQE QLQKEQEKAF FAQFQLDEET GEFLPIQPAQ HIQTDTSGSA SYSQVAHIPK QDALYFEDCM QLLAETFPFV DDHESLALDI PSHAESSVFT APHQAQSLNS SLEAAMTDLS SIEQDMEQVW QELFSIPELQ CLNTENKQLA DTTAVPSPEA TLTEMDSNYH FYSSISSLEK EVGNCGPHFL HGFEDSFSSI LSTDDASQLT SLDSNPTLNT DFGDEFYSAF IAEPDGGSM PSSAAISQSL SELLDGTIEG CDLSLCKAFN PKHAEGTMEF NDSDSGISLN TSPSRASPEH SVESSIYGDP PPGFSDSEME ELDSAPGSVK QNGPKAQPAH SPGDTVQPLS PAQGHSAPMR ESQCENTTKK EVPVSPGHQK APFTKDKHSS RLEAHLTRDE LRAKALHIPF PVEKIINLPV DDFNEMMSKE QFNEAQLALI RDIRRRGKNK VAAQNCRKRK LENIVELEQD LGHLKDEREK LLREKGENDR NLHLLKRRLS TLYLEVFSML RDEDGKPYSP SEYSLQQTRD GNVFLVPKSK KPDTKKN
<b>Source</b>	Baculovirus
<b>Target Names</b>	Nfe2l2
<b>Protein Names</b>	Recommended name: Nuclear factor erythroid 2-related factor 2 Short name= NF-E2-related factor 2 Short name= NFE2-related factor 2 Alternative name(s): Nuclear factor, erythroid derived 2, like 2
<b>Expression Region</b>	1-597
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.