



Recombinant Mouse Myocyte-specific enhancer factor 2A (Mef2a)

Product Code	CSB-YP726688MO
Abbreviation	Mef2a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60929
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MGRKKIQITR IMDERNRQVT FTKRKFGMLK KAYELSVLCD CEIALIIFNS SNKLFQYAST DMDKVLKYT EYNEPHESRT NSDIVETLRK KGLNGCESPD ADDYFEHSPL SEDRFSKLE DSDFIFKRG PGLPPQNFMS SVTVPVTSPN ALSYTNPGSS LVSPSLAASS TLADSSMLSP PPATLHRNVS PGAPQRPPST GSASGMLSTT DLTVPNGAGN SPVGNFVNS RASPNIIGHT GANSLGKVMP TKSPPPGGG SLGMNSRKP LRVVIPPSSK GMMPLSEEE EELNAQRIS SSQATQPLAT PVVSVTTPSL PPQGLVYSAM PTAYNTDYSL TSADLSALQG FTSPGMLSLG QASAWQQHHL GQAALSSLVA GGQLSQGSNL SINTNQNINI KSEPISPPRD RMTPSGFQQQ QQQPQQQPPP QPPQPQRQE MGRSPVDSLS SSSSSYDGSD REDPRGDFHS PIVLGRPPNT EDRESPSVKR MRMDTWVT
Source	Yeast
Target Names	Mef2a
Protein Names	Recommended name: Myocyte-specific enhancer factor 2A
Expression Region	1-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.