



Recombinant Methylococcus capsulatus Malate dehydrogenase (mdh)

Product Code	CSB-BP726702MFM
Abbreviation	mdh
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60B71
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MKTPVHVAVT GAAGQIAYSL LFRIAVGDLF GPHQPVLKL LDVPSAERVL EGVAMELDDC ASPLLQEIEV SSDPAEVFDG AEAVFMLGAT PRGPGMERRD LLQVNADIFS AQGRALNESA SRRVKILVVG NPANTNALIA QRNAPDLAPG CFSAMTRLDH NRATSLLARH CGCNVAEISR VVIWGNHSPT QYPDLHHARV KGKPALSLVD PAWYVETFIP TVQQRGASVI AIRGKSSAAS AANAALDHMR SWFLGTPKDD WVSMTVSSDG SYGIAEGLMF SFPVTIENGR FRIVQDLPLD TFSRERLRLT EVELLEERAM VSHLL
Source	Baculovirus
Target Names	mdh
Protein Names	Recommended name: Malate dehydrogenase EC= 1.1.1.37
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.