



Recombinant *Methylococcus capsulatus* Protein RecA (recA)

Product Code	CSB-YP726707MFM
Abbreviation	recA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60BS8
Product Type	Recombinant Protein
Immunogen Species	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MDENKRKALG AALSQIEKQF GKGSVMRMGD VAAVRDIDVI PTGSLALDIA LGCGGLPRGR IVEIYGPESS GKTTLTLEVI AQAQKTGGVA AFVDAEHALD PVYAEKIGVN LDDLLVSQPD TGEQALEIAD VLVRS GGVDV VVIDSVAALT PKAELEGEMG DSHMGLQARL MSQALRKLTA NIKRSNTLVI FINQIRMKIG VMFGNPETTT GGNALKFYAS VRLDIRRSGA IKNGDEVIGN ETRVKVVKNK VAPPFRTADF EILYGEGISR EGELIDLGVN HDIVQKSGSW YSYGGDRIGQ GKDNVRIYLK EHPEVAAAIE AAVREKALAG FQHAPSAERA ALEESGF
Source	Yeast
Target Names	recA
Protein Names	Recommended name: Protein RecA Alternative name(s): Recombinase A
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.