



Recombinant Methylococcus capsulatus 3-dehydroquinate synthase (aroB)

Product Code	CSB-MP726709MFM
Abbreviation	aroB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60BY2
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MKTLHVELGE RGYPIYIGRG LLGHPDLIQA HLPGGQVLVV TNEVVAPLYL DRMLASLAGK DTGSVVLDPG EAHKTLD SAM AVFDALLARR FGRNAAIVAL GGGVIGDLAG FAAACYQRGV PFIQVPTTLL SQVDSSVGGK TAVNHPRGKN MIGAFYQPRC VLADTDLTLD LPDRELSAGL AEVIKYGFIR DPEFLAWLEA NVERLLQRDP EALAYAIERS CINKAEIVAE DETETGVRAT LNLGHTFGHA IETGMGYGVC LHGEAVAIGM CQAADLSRRL GWIGDDEVAR VIRLLERARL PVVPPRELDA DAFLEHMAVD KKNVDGGLRL VLLKSLGEAT LPVAVDAGLL RATLECYGR
Source	Mammalian cell
Target Names	aroB
Protein Names	Recommended name: 3-dehydroquinate synthase EC= 4.2.3.4
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.