



# Recombinant Cat Growth arrest and DNA damage-inducible protein GADD45 alpha (GADD45A)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP726719CA  |
| <b>Abbreviation</b>      | GADD45A   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q60G15  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Felis catus (Cat) (Felis silvestris catus)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MTLEEFSAAGE QKTERMDKVG DALEEVLSKA LSQRTITVGV YEAAKLLNVD<br>PDNVVLCLLA ADEDDDRDVA LQIHFTLIQA FCCENDINIL RVSNPGRLEA<br>LLLLETDASP AASEGAEQPP DLHCVLVTNP HSSQWKDPAL SQLICFCRES<br>RYMDQWVPVI NLPER   |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | GADD45A   |
| <b>Protein Names</b>     | Recommended name: Growth arrest and DNA damage-inducible protein<br>GADD45 alpha  |
| <b>Expression Region</b> | 1-165   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Target Details</b>    | This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. This protein responds to environmental stresses by mediating activation of the p38/JNK pathway via MTK1/MEKK4 kinase. The DNA damage-induced transcription of this gene is mediated by both p53-dependent and -independent mechanisms. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.