



Recombinant Mouse NGFI-A-binding protein 1 (Nab1)

Product Code	CSB-MP726736MO
Abbreviation	Nab1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61122
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MATALPRTLGL ELQLYRILQK ANLLSYFDAF IQQGGDDVQQ LCEAGEEEFL EIMALVGMAK KPLHVRRLQK ALRDWVTNPG LFNQPLTSLP VSSIPIYKLP EGSPTWLGIS CNSYERSSSS REPHLKIPKC AATTCVQSLG QGKSEVGSLA LQSVSDSRLW QGHHATESEH SLSPADLGSP ASPKESSEAL DAAAALSVAE CVERMAPTLP KSDLSEVKEL LKNNKKLAKM IGHIFEMSDE DPHKEEEIRK YSAIYGRFDS KRKDGKHLTL HELTVNEAAA QLCVKDNALL TRRDELFALA RQVSREVTYK YTYRTRLKLC GERDELSPKR IKIEDGFPDF QESVPTLFQQ ARAKSEELAG LGSQQAEEKGM AKQMELLCAQ AGYERLQQR RLTAGLYRQS SGEQSPDGGL PSDSSDGQGE RPLNLRIPSV QNRQPHHFVV DGELSRLYSS EAKSHSSESL GILKDYPHSA FTLEKKVIKT EPEDSR
Source	Mammalian cell
Target Names	Nab1
Protein Names	Recommended name: NGFI-A-binding protein 1 Alternative name(s): EGR-1-binding protein 1
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.