



Recombinant Mouse T-complex protein 1 subunit zeta-2 (Cct6b)

Product Code	CSB-EP726745MO-B
Abbreviation	Cct6b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61390
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAAIKIANPG AEVTRSQAAL AVNICAARGL QDVLRLPTLGP KGALKMLVSG AGDIKLT KD G NVLLHEMQIQ HPTASIIAKV AAAQDHVTGD GTTSNVLIIG ELLKQADLYI SEGLHPRIIT EGFDVAKTKA LEVLDEIKVQ KEMKREILLD VARTSLQTKV HAELADILTE AVVDSVLAIR RPGVPIDLFM VEIVEMRHKS ETDTQLIRGL VLDHGARHPR MRKQVRDAYI LTCNVSLEYE KTEVSSGFFY KTVEEKEKLV KAERKFIEDR VQKIIDLKQK VCAESNKG FV VINQKGIDPV SLEMLAKHNI VALRRAKRRN LERLTLACGG LAVNSFEGLS EECLGHAGLV FEYALGEEKF TFIEDCVNPL SVTLLVKGP N KHTLIQIKDA LRDGLRAVKN AIEDGCVVPG AGAVEVAIAE ALVNYKHRVQ GRVRLGIQAF ADALLIIPKV LAQNSGYDLQ ETLIKIQT KH AESKELLGID LNTGEPMAAA EAGIWDNYCV KKHLLHSCTV IATNILLVDE IMRAGMSSLR D
Source	E.coli
Target Names	Cct6b
Protein Names	Recommended name: T-complex protein 1 subunit zeta-2 Short name= TCP-1-zeta-2 Alternative name(s): CCT-zeta-2 Short name= Cctz-2
Expression Region	1-531
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.