



Recombinant Mouse Fascin (Fscn1)

Product Code	CSB-MP726751MO
Abbreviation	Fscn1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q61553
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	TANGTAEAV QIQFGLISCG NKYLTAEAFG FKVNASASSL KKKQIWTLEQ PPDEAGSAAV CLRSHLGRYL AADKDGNVTC EREVPDGDGR FLVVAHDDGR WSLQSEAHRR YFGGTEDRLS CFAQSVSPAE KWSVHIAMHP QVNIYSVTRK RYAHLSPARPA DEIAVDRDVP WGVDSLITLA FQDQRYSVQT SDHRFLRHDG RLVARPEPAT GFTLEFRSGK VAFRDCEGRY LAPSGPSGTL KAGKATKVGK DELFALEQSC AQVVLQAANE RNVSTRQGMD LSAHQDEETD QETFQLEIDR DTRKCAFRTTH TGKYWTLTAT GGVQSTASTK NASCYFDIEW CDRRITLRAS NGKFVTAKKN GQLAASVETA GDSEFLMKL INRPIIVFRG EHGFIGCRKV TGTLDANRSS YDVFQLEFND GAYNIKDSTG KYWTVGSDSS VTSSSDTPVD FFLEFCDYNK VALKVGGRYL KGDHAGVLKA CAETIDPASL WEY
Source	Mammalian cell
Target Names	Fscn1
Protein Names	Recommended name: Fascin Alternative name(s): Singed-like protein
Expression Region	2-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>